Applied Thermal Control Ltd

39 Hayhill Industrial Estate Barrow upon Soar, Leicestershire LE12 8LD, United Kingdom

UPDATED GUIDANCE USE OF TAP WATER IN ATC CHILLERS



Date 12/06/2024 Page 1

Earlier documentation may have indicated that the use of tap water in ATC chillers was permissible under certain conditions. Based on further evaluation and our ongoing commitment to system reliability, we now advise against the use of tap water in any ATC chiller models. This updated recommendation reflects a more comprehensive understanding of global water quality variation.

WHY THIS MATTERS

Water quality can vary widely by region and application. Tap water may contain:

- High mineral content (e.g. calcium, magnesium)
- pH imbalance
- Biological matter (e.g. algae, bacteria)
- Chemical additives (e.g. chlorine, fluoride)

These factors can contribute to:

- Internal corrosion
- Degradation of seals and materials
- Biological fouling
- Reduced heat transfer performance

As such, continued use of unapproved fluids, including untreated tap water, may affect long-term performance in some instances.

CURRENT RECOMMENDATION

To support the ongoing reliability and efficiency of ATC chillers, we recommend using only ATC-approved fluids.

PREFFERED FLUID

Hexid - Providing corrosion protection, freeze protection, algae inhibition, and effective thermal performance.

If Hexid is not available in your region, please contact ATC for guidance.


SUPPORT

For questions regarding coolant compatibility or fluid selection for your specific system, please contact ATC Technical Support.

Installation, Operation & Service Manual **G04/K04**

Rev 1 - 05/04/23





Installation, Operation & Service Manual **G04/K04**

Rev 1 - 05/04/23

Table of Contents

1.0	INTRODUCTION	3
1.1	SAFETY NOTICES	4
1.2	UNPACKING UNITS WEIGHING OVER 60kg(133lbs) ON CASTORS	5
1.3	SITE REQUIREMENTS	5
1.4	WARRANTY REGISTRATION	7
2.0	PRODUCT SPECIFICATION	8
3.0	QUICK START GUIDE	9
4.0	INSTALLATION	10
4.1	HOSE RECOMMENDATIONS	10
4.2	CONNECTING ADAPTERS TO PRODUCT BULKHEAD FITTINGS	10
4.3	BACKFILLING	11
4.4	VOLTAGE SELECTION	11
4.4	FILLING PROCEDURES	11
4.5	DRAINING A SYSTEM	12
5.0	OPERATION	13
5.1	CHANGING SET POINT	14
5.2	REGULATING PRESSURE/FLOW	14
	ASIC TROUBLESHOOTING	
7.0 EN	ND USER MAINTANANCE	17
8.0 DE	ECLARATIONS AND APPROVAL	18

Installation, Operation & Service Manual **G04/K04**

Rev 1 - 05/04/23

1.0 INTRODUCTION

By selecting a K04/G04 series chiller you have invested in many years of experience in the design and manufacture of precision temperature control instrumentation.

ATC has built your chiller without compromise to meet the objectives of performance and reliability. Please read this manual carefully to ensure you understand the operation of the machine and how to use the unit safely and efficiently.

If you have any questions regarding installation or repair of this unit, please contact ATC direct.

Applied Thermal Control Ltd.

39 Hayhill Industrial Estate, Barrow Upon Soar

Leicestershire,

LE12 8LD

United Kingdom

Tel: +44 (0) 1530 839998

E-mail: support@app-therm.com

Installation, Operation & Service Manual **G04/K04**

Rev 1 - 05/04/23

1.1 SAFETY NOTICES

GENERAL SAFETY INFORMATION

For your safety, we draw your attention to the following warning and caution marks throughout the manual. Warning symbols can be found on the unit. Ensure you have read through all warnings before starting the unit.

The safe operation of ATC products always remains the responsibility of the operator. This equipment is intended to be used as a liquid temperature conditioning device – it requires no external pump, nor any further manipulation of temperature. If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired. Failure to comply with a 'warning' may result in personal injury or death. ATC does not accept any liability for injury caused through use of this equipment.

-	d. Failure to comply with a 'warning' may result in personal injury or death. ATC does ept any liability for injury caused through use of this equipment.
0	Caution: Failure to comply with a caution will invalidate product warranty and absolve ATC from any liability, howsoever caused, and could result in permanent damage to equipment.
0	Caution: Filling/topping up of the tank should only be undertaken with the unit switched off, to prevent back-filling of the fluid.
0	Caution: This product contains no user-serviceable parts. Repair and service requires specialized knowledge and tools to be provided by ATC or its local agent. Any unauthorized tampering with the heat exchanger system automatically invalidates warranty.
lack	Warning: Hot and cold surfaces are present during operation. Take caution and care when touching pump during operation.
A	Warning: Water pressures of up to 10 bar during operation.
A	Warning; Water and electricity near one another. Always ensure the unit is isolated before service. The product is protected from overcurrent by fuses. Never bypass this component.
A	During fault diagnostics and maintenance, it may be necessary to remove panels, which expose the operator to the dangers of pressurized systems, hot or cold pipes and electrical circuits. Only qualified personnel who are aware and equipped to deal with these systems should only carry out such work.
lack	Any temporary electrical supply to the chiller should be correctly earthed and connected through an earth leakage trip.
A	In case of unexpected coolant leakage, safety glasses should always be worn when the chiller is operated with the covers removed.
Λ	Under no circumstances leave the cooler unattended with the side panels removed.

Never alter settings of pressure switches, overloads, circuit breakers or any safety

device without first consulting Applied Thermal Control.

Installation, Operation & Service Manual **G04/K04**

Rev 1 - 05/04/23

1.2 UNPACKING UNITS WEIGHING OVER 60kg(133lbs) ON CASTORS

Please check that both the packaging and the unit are undamaged. If there is any doubt, it is vital that you inform both ATC and the carrier. There are no hidden shipping bolts or other fixings. You should inspect the packaging for signs of transit damage before signing for the unit, and if possible, unpack the unit before signing. Once you have signed for the goods, ATC cannot be held responsible for any transit damage subsequently found.

As the unit weighs >60kg, ATC recommends it should be lifted with slings through the underside, using a forklift or overhead crane. ATC highly recommends that it is not manually lifted, and that safe slinging and lifting practices are adhered to.

Remove the unit from its original packaging and ensure that there is no packaging left around the cooling ducts. There is no internal product packaging that requires the chiller to be opened.

Please retain all packaging in the unlikely event that the chiller needs to be returned to our local representatives.

1.3 SITE REQUIREMENTS

5%. Before starting storing outside <i>operating</i> C to +65°C. 1°C (+88°F), decreasing ature. filling and allowing air to
Storing outside <i>operating</i> C to +65°C. 1°C (+88°F), decreasing ature.
C to +65°C. 1°C (+88°F), decreasing ature.
1°C (+88°F), decreasing ature.
1°C (+88°F), decreasing ature.
ature.
filling and allowing air to
it breaker is rated at 15
circuit breaker is also
se overload handle in
g recommended bending
ved.
, cleanest air that can be
s air in from the left-
ect heated air from the
h the fluid to be used.
er-glycol mixtures such as
dling the nominal



Installation, Operation & Service Manual **G04/K04**

Rev 1 - 05/04/23

9	Indoor use only. Altitude up to 2000m above sea level. Ensure the unit has adequate						
	ventilation.						
	Caution: Always use ATC recommended fluids in our products – many sealing						
J	compounds and materials are present and unapproved fluids have the potential to						
	corrode your application and damage seals.						
•	Caution: Do not use inadequately rated wiring. Consult an electrician if you are unsure.						
	Caution: The safety of any system incorporating the equipment is the responsibility of						
J	the assembler of the system.						
	Caution; Do not replace detachable mains cords with inadequately rated cords. Contact						
	ATC for appropriately rated products.						

Installation, Operation & Service Manual **G04/K04**

Rev 1 - 05/04/23

1.4 WARRANTY REGISTRATION

Please visit the website warranty registration page to ensure ATC can offer you the best possible support.

https://www.app-therm.com/warranty-registration/

a) For how long is my ATC product under warranty?

ATC provides a comprehensive return to base 2-year parts, 1-year labour warranty from delivery as standard on all new equipment, provided it has been installed and operated in accordance with the manual.

b) Where will ATC fulfill the product warranty?

ATC's standard warranty terms are Return to Base (RTB) – issues with chillers are often easily solvable over the phone or email, or by reviewing ATC's technical guidance on the web and in the product manual. On occasion, at the discretion of ATC, goods may be serviced on site FOC or a service loan unit may be supplied. Warranty cover excludes the cost of travel by engineers and loan unit rental charges. Obtaining onsite service for a product, even in full warranty, is a chargeable service.

- Who is liable for shipping charges in the event of warranty failure?
 During the first year of the warranty period, freight costs to and from ATC are covered by ATC.
 During the second year of the warranty, freight costs to and from ATC are payable by the customer.
- d) I'm experiencing problems with my chiller. It's within warranty what do I do next? Contact ATC to discuss the issue you are having on +44(0)1530 839998 or support@app-therm.com. Be sure to have your model number and serial number on-hand to aid those attempting to solve remotely.

e) Telephone support couldn't fix my chiller – what do I do next?

An RMA form must be completed. This allows both the end-user and ATC to clarify your details, to set the party responsible for shipping costs, and to set a different return address if desired. Shipping advice is provided, and the end-user must sign a declaration that states the unit is safe to handle. Return the form by email for fastest response.

- f) What happens if my chiller failed outside warranty or requires non-warranty repair work? A purchase order will be requested to cover an initial inspection this will only be invoiced if the inspection shows there is no fault. If packaging is required, i.e. a crate, a separate charge will be levied. If the end user prefers ATC to arrange a collection, a shipping charge may be levied.
- g) Our process must continue running can we have a loan unit whilst our chiller is in repair?

 ATC hold several standard air-cooled chillers at the factory for the sole purpose of offering for loan.

 These are available on a first come, first-serve basis. Models up-to 3kW capacity are available.

Installation, Operation & Service Manual **G04/K04**

Rev 1 - 05/04/23

2.0 PRODUCT SPECIFICATION

Weight 198lbs 1600 x W513 x H888mm 123.6 x W20.2 x H35" 336lk @1m 1000 10	Model	G04/K04
Dimensions L600 x W513 x H888mm L23.6 x W20.2 x H35" Sold So	Weight	<u> </u>
Noise level Toolless access Technology Vapor Compression ShW (Ambient temperature at +30°C with setpoint of +20°C) Refrigerant Refrigerant Control method Control method PID via Hot Gas Bypass Stability Resolution Sensor type System volume PD pump options Pressure relief Process fittings Process fittings Supply (-0spec) Supply (-0spec) Supply (-12 pec) Temp range alarm Low fluid flow alarm Low fluid flow alarm Temp range alarm Fridge HP switch Motor thermal cutout Overcurrent protection Overcurrent restart Emergency off Not present Supon Vapor Compression Shabilty 1,942BTU/h 0,995TR G04: R290 (Propane) K04: R407C Hove Too Hot Aud Too Possor (extended) PID via Hot Gas Bypass PID via Hot Gas Bypass PID via Hot Gas Bypass PTIDO Sensor type PT100 System volume 3L PD pump options PF, P10, P17 Internal, PRV, 20-150psi 1/2" BSPPF Supply (-0spec) 230Vac 50Hz 1~9A Supply 1-2spec) 220Vac 60Hz 1~/2~ 9A Supply (-12 pec) Standard Standard Standard Standard Standard Standard Standard Standard, via MCB Not present 2-year parts, 1-year labor Enhanced on request	Form factor	Floorstanding
Toolless access Technology Cooling capacity (Ambient temperature at +30°C with setpoint of +20°C) Refrigerant Refrigerant Temperature range Control method Control method Sensor type System volume PD pump options Pressure relief Process fittings Supply (-0spec) Supply (-2spec) Supply (-2spec) Supply (-1spec) Temp range alarm Fridge HP switch Motor thermal cutout Overcurrent protection Overcurrent restart Emergency off Warranty options No Vapor Compression 3.5kW 11,942BTU/h 0.995TR 604: R290 (Propane) K04: R40°C ENGLE R40°C 11,942BTU/h 10,942BTU/h 10,94BTU/h 10,942BTU/h 10,	Dimensions	
Toolless access Technology Cooling capacity (Ambient temperature at +30°C with setpoint of +20°C) Refrigerant Refrigerant Temperature range -20°C to +65°C (extended) Control method Control method Sensor type System volume PD pump options Process fittings Process fittings Supply (-0spec) Supply (-2spec) Supply (-2spec) Temp range alarm Low fluid flow alarm Temp range alarm Fridge HP switch Motor thermal cutout Overcurrent protection Overcurrent restart Emergency off Warranty options No SabW 11,942BTU/h 0.995TR 604: R290 (Propane) K04: R407C 11,942BTU/h 0.995TR 11,942BTU/h 0.995TR 2004 R290 (Propane) K04: R407C 11,942BTU/h 0.995TR 2004: R290 (Propane) K04: R407C 11,942BTU/h 0.995TR 2004: R407C 11,942BTU/h 0.995TR 2004: R407C 11,942BTU/h 0.995TR 2004: R290 (Propane) K04: R407C FNO.995TR 2004: R407C 11,942BTU/h 0.995TR 2004: R290 (Propane) K04: R407C 11,942BTU/h 0.995TR 20,940 R04: R290 (Propane) K04: R407C FNO.995TR 604: R290 (Propane) K04: R407C FNO.995TR 604: R290 (Propane) K04: R407C FNO.995TR 604: R290 (Propane) K04: R407C FNO.995TR 604: R290 (Propane) K04: R407C FNO.996TR 11,942BTU/h 0.995TR 604: R407C FNO.996TR 604: R40C FNO.996TR 604: R40C FNO.996TR 60	Noise level	53dB(A)@1m
Technology Cooling capacity (Ambient temperature at +30°C with setpoint of +20°C) Refrigerant Refrigerant Control method Control method Stability Resolution Sensor type PD pump options Process fittings Process fittings Supply (-0spec) Supply (-2spec) Supply (-2spec) Low fluid flow alarm Temp range alarm Fridge HP switch Motor thermal cutout Overcurrent protection Overcurrent protection Warranty options Refrigerant A.SkW 11,942BTU/h 0.995TR 604: R290 (Propane) K04: R407C +4°C to +35°C (standard) -20°C to +65°C (extended) PID via Hot Gas Bypass ±0.1°C PT100 31. PP PT100 31. PP PT100 31. PP P, P10, P17 Internal, PRV, 20-150psi 1/2" BSPPF 230Vac 50Hz 1~9A 220Vac 60Hz 1~/2~ 9A CE, UKCA Visual, lamp Visual, lamp Visual, lamp Temp range alarm Fridge HP switch Motor thermal cutout Overcurrent protection Overcurrent restart Emergency off Warranty options		
Cooling capacity (Ambient temperature at +30°C with setpoint of +20°C) Refrigerant Resolution PID via Hot Gas Bypass Stability ±0.1°C Resolution O.1°C Sensor type PT100 System volume PD pump options Pressure relief Pressure relief Process fittings Process fittings Process fittings Supply (-0spec) Supply (-0spec) Supply (-2spec) PT100 Supply (-2spec) Supply (-2spec) Supply (-2spec) Supply (-2spec) Supply (-3pec) Supply	Technology	Vapor Compression
K04: R407C		11,942BTU/h
-20°C to +65°C (extended) Control method BID via Hot Gas Bypass Stability ±0.1°C Resolution Sensor type PT100 System volume BPD pump options PF, P10, P17 Pressure relief Process fittings Process fittings Process fittings BSPPF Supply (-0spec) Supply (-2spec) Supply (-2spec) Brote Andrew Supply (-2spec) CE, UKCA Empty fluid tank alarm Emprange alarm Visual, lamp Temp range alarm Fridge HP switch Motor thermal cutout Overcurrent protection Overcurrent restart Emergency off Warranty options PTC PD pump option BL PT00 BSPPF Supply (-2spec) BSPPF BSUpply (-2spec) BSPPF BSPPF BSPPF BSPPF BSUpply (-2spec) BSPPF BSPPF BSPPF BSPPF BSUpply (-2spec) BSPPF BSPPF BSPPF BSUpply (-2spec) BSPPF BSPPF BSUpply (-2spec) BSPF BSPF BSPF BSPF BSPF BSPF BSPF BSPF	Refrigerant	K04: R407C
Control method Stability ±0.1°C Resolution Sensor type PT100 System volume BPD pump options PPD pump options Process fittings	Temperature range	+4°C to +35°C (standard)
Stability ±0.1°C Resolution 0.1°C Sensor type PT100 System volume 3L PD pump options P5, P10, P17 Pressure relief Internal, PRV, 20-150psi Process fittings 1/2" BSPPF Supply (-0spec) 230Vac 50Hz 1~9A Supply (-2spec) 220Vac 60Hz 1~/2~ 9A 1st party approvals CE, UKCA Empty fluid tank alarm Visual, lamp Low fluid flow alarm Visual, lamp Temp range alarm Visual, lamp Fridge HP switch Standard Motor thermal cutout Standard Overcurrent protection Overcurrent restart Manual Emergency off Not present Warranty options Versure Resolution Outlied 10.1°C Solution United Standard Standard Not present 2-year parts, 1-year labor Enhanced on request		-20°C to +65°C (extended)
Resolution 0.1°C Sensor type PT100 System volume 3L PD pump options P5, P10, P17 Pressure relief Internal, PRV, 20-150psi Process fittings 1/2" BSPPF Supply (-0spec) 230Vac 50Hz 1~ 9A Supply (-2spec) 220Vac 60Hz 1~/2~ 9A 1st party approvals CE, UKCA Empty fluid tank alarm Visual, lamp Low fluid flow alarm Visual, lamp Temp range alarm Visual, lamp Fridge HP switch Standard Motor thermal cutout Standard Overcurrent protection Standard, via MCB Overcurrent restart Manual Emergency off Not present Warranty options 2-year parts, 1-year labor Enhanced on request	Control method	
Sensor type PT100 System volume 3L PD pump options P5, P10, P17 Pressure relief Internal, PRV, 20-150psi Process fittings 1/2" BSPPF Supply (-0spec) 230Vac 50Hz 1~ 9A Supply (-2spec) 220Vac 60Hz 1~/2~ 9A 1st party approvals CE, UKCA Empty fluid tank alarm Visual, lamp Low fluid flow alarm Visual, lamp Temp range alarm Visual, lamp Fridge HP switch Standard Motor thermal cutout Standard Overcurrent protection Standard, via MCB Overcurrent restart Manual Emergency off Not present Warranty options 2-year parts, 1-year labor Enhanced on request	Stability	
System volume PD pump options PF, P10, P17 Pressure relief Internal, PRV, 20-150psi Process fittings Process fittings Process fittings Supply (-0spec) Supply (-2spec) 220Vac 50Hz 1~9A Supply (-2spec) 220Vac 60Hz 1~/2~9A CE, UKCA Empty fluid tank alarm Visual, lamp Low fluid flow alarm Visual, lamp Temp range alarm Visual, lamp Fridge HP switch Standard Motor thermal cutout Standard Overcurrent protection Overcurrent restart Emergency off Warranty options System volume System volume 1/2" BSPPF Supply (-2spec) 220Vac 50Hz 1~9A Supply (-2spec) Standar 1~/2~9A Supply (-2spec) Supply (-2	Resolution	0.1°C
PD pump options P5, P10, P17 Pressure relief Internal, PRV, 20-150psi Process fittings 1/2" BSPPF Supply (-0spec) 230Vac 50Hz 1~ 9A Supply (-2spec) 220Vac 60Hz 1~/2~ 9A 1st party approvals CE, UKCA Empty fluid tank alarm Visual, lamp Low fluid flow alarm Visual, lamp Temp range alarm Visual, lamp Fridge HP switch Standard Motor thermal cutout Standard Overcurrent protection Standard, via MCB Overcurrent restart Manual Emergency off Not present Warranty options P5, P10, P17 Internal, PRV, 20-150psi 1/2" BSPPF Supply (-0spec) 230Vac 50Hz 1~ 9A CE, UKCA Emuty Standard Visual, lamp Fridge HP switch Standard Standard Overcurrent protection Standard Overcurrent protection Standard, via MCB Overcurrent restart Manual Emergency off Not present 2-year parts, 1-year labor Enhanced on request	Sensor type	PT100
Pressure relief Internal, PRV, 20-150psi Process fittings 1/2" BSPPF Supply (-0spec) 230Vac 50Hz 1~ 9A Supply (-2spec) 220Vac 60Hz 1~/2~ 9A 1st party approvals CE, UKCA Empty fluid tank alarm Visual, lamp Low fluid flow alarm Visual, lamp Temp range alarm Visual, lamp Fridge HP switch Standard Motor thermal cutout Standard Overcurrent protection Standard, via MCB Overcurrent restart Manual Emergency off Not present Warranty options University 1/2" BSPPF 2004 2004 2004 2004 2004 2004 2004 200	System volume	
Process fittings 1/2" BSPPF Supply (-0spec) 230Vac 50Hz 1~ 9A Supply (-2spec) 220Vac 60Hz 1~/2~ 9A 1st party approvals CE, UKCA Empty fluid tank alarm Visual, lamp Low fluid flow alarm Visual, lamp Temp range alarm Visual, lamp Fridge HP switch Standard Motor thermal cutout Standard Overcurrent protection Standard, via MCB Overcurrent restart Manual Emergency off Not present Warranty options Warranty options		P5, P10, P17
Supply (-0spec) 230Vac 50Hz 1~ 9A Supply (-2spec) 220Vac 60Hz 1~/2~ 9A 1st party approvals CE, UKCA Empty fluid tank alarm Visual, lamp Low fluid flow alarm Visual, lamp Temp range alarm Visual, lamp Fridge HP switch Standard Motor thermal cutout Standard Overcurrent protection Standard, via MCB Overcurrent restart Manual Emergency off Not present Warranty options 230Vac 50Hz 1~ 9A 200Vac 60Hz 1~/2~ 9A CE, UKCA Visual, lamp Visual, lamp Standard Standard Notor thermal cutout Standard Overcurrent protection Standard Overcurrent restart Manual Emergency off Not present 2-year parts, 1-year labor Enhanced on request	Pressure relief	•
Supply (-2spec) 220Vac 60Hz 1~/2~ 9A 1st party approvals CE, UKCA Empty fluid tank alarm Visual, lamp Low fluid flow alarm Visual, lamp Temp range alarm Visual, lamp Fridge HP switch Standard Motor thermal cutout Standard Overcurrent protection Standard, via MCB Overcurrent restart Manual Emergency off Not present Warranty options 2-year parts, 1-year labor Enhanced on request	Process fittings	-
1st party approvals CE, UKCA Empty fluid tank alarm Visual, lamp Low fluid flow alarm Visual, lamp Temp range alarm Visual, lamp Fridge HP switch Standard Motor thermal cutout Standard Overcurrent protection Standard, via MCB Overcurrent restart Manual Emergency off Not present Warranty options CE, UKCA CHANCE Visual, lamp Standard Visual, lamp Standard Standard Standard Overcurrent protection Standard Overcurrent restart Manual Emergency off Not present 2-year parts, 1-year labor Enhanced on request		230Vac 50Hz 1~ 9A
Empty fluid tank alarm Low fluid flow alarm Temp range alarm Fridge HP switch Motor thermal cutout Overcurrent protection Overcurrent restart Emergency off Warranty options Visual, lamp Visual, lamp Standard Visual, lamp Standard Standard Motor thermal cutout Standard, via MCB Manual Emergency off Vot present 2-year parts, 1-year labor Enhanced on request		220Vac 60Hz 1~/2~ 9A
Low fluid flow alarm Temp range alarm Visual, lamp Fridge HP switch Standard Motor thermal cutout Standard Overcurrent protection Overcurrent restart Emergency off Warranty options Visual, lamp Standard Standard MCB Overcurrent restart Manual Emergency off Vot present 2-year parts, 1-year labor Enhanced on request	1 st party approvals	CE, UKCA
Temp range alarm Visual, lamp Fridge HP switch Standard Motor thermal cutout Standard Overcurrent protection Standard, via MCB Overcurrent restart Manual Emergency off Not present 2-year parts, 1-year labor Enhanced on request		Visual, lamp
Fridge HP switch Standard Motor thermal cutout Standard Overcurrent protection Standard, via MCB Overcurrent restart Manual Emergency off Not present Warranty options 2-year parts, 1-year labor Enhanced on request	Low fluid flow alarm	Visual, lamp
Motor thermal cutout Standard Overcurrent protection Standard, via MCB Overcurrent restart Manual Emergency off Not present 2-year parts, 1-year labor Enhanced on request	Temp range alarm	Visual, lamp
Overcurrent protection Standard, via MCB Overcurrent restart Manual Emergency off Not present 2-year parts, 1-year labor Enhanced on request	Fridge HP switch	Standard
Overcurrent restart Manual Emergency off Not present 2-year parts, 1-year labor Enhanced on request	Motor thermal cutout	Standard
Emergency off Not present 2-year parts, 1-year labor Enhanced on request	Overcurrent protection	Standard, via MCB
Warranty options 2-year parts, 1-year labor Enhanced on request	Overcurrent restart	Manual
Warranty options Enhanced on request	Emergency off	Not present
	Warranty options	
Rated duty cycle SI Continuous	Rated duty cycle	S1 Continuous



Installation, Operation & Service Manual **G04/K04**

Rev 1 - 05/04/23

3.0 QUICK START GUIDE

	QUICK START GUIDE	
1	Remove all the packaging (Section 1.2)	
2	Connect chiller's outlet and inlet process to your application	
3	Fill up the chiller with appropriate coolant	
4	Use the name plate on the back of the chiller to verify the appropriate	
	voltage and supply the power using provided cables.	
5	Turn the power switch ON	
6	Slowly add more coolant as process lines fill and the coolant is	
	circulated around the process system.	
7	Using the controller, set the desired temperature.	

Installation, Operation & Service Manual **G04/K04**

Rev 1 - 05/04/23

4.0 INSTALLATION

4.1 HOSE RECOMMENDATIONS

Having ensured that your installation meets all site requirements, it is best practice that the fluid lines between your application and the chiller have the following characteristics.

- **Short in length** this reduces friction-based pressure drop and addition ambient heat load.
- 2 Large diameter bore at least 12mm (1/2").
- Free from 90° bends to limit the effects of water hammer. If this cannot be avoided, sharp changes of direction should be minimized so far as possible. Doing this correctly can yield higher pump performance and extend time between maintenance intervals. It will also reduce electrical energy consumption.
- Clean If your installation is to existing pipe work, it is good practice to flush the system with either a commercially available central heating cleaner or 5% acetic acid solution. The system should be flushed clean with tap water to remove all traces of cleaner prior to filling the system. Failing this, it is recommended to use a domestic bleach in solution with tap water, diluted to the point where the bleach can longer be smelled by human nose.
- **Opaque, ideally black** to inhibit light passing through the tube and algae building up. Alternatively, solid ABS or copper pipe can be used where application chemistry allows.
- Insulation, where low temperature process is planned the process line from chiller to application contains the feed of low temperature fluid. Insulation prevents heat from entering this line and can promote better stability. Uninsulated return lines are helpful where free cooling can be obtained by allowing heat to transfer to air likewise, insulating the return line is helpful if the fluid temperature is below ambient.



Caution; Never use transparent tubing. UV light will pass through, prompting growth of organic contamination.

4.2 CONNECTING ADAPTERS TO PRODUCT BULKHEAD FITTINGS

- Standard units are supplied with 1/2" British Standard Pipe Parallel Female (BSPPF) threads (also known as G threads (ISO228)) by default. These fittings are not valved and will 'drop' the volume of the system if left open to atmosphere.
- 2 Ensure the appropriate thread sealants are used in the fitting of adapters to hose. For metallic mating parts, we recommend Loctite 577. For plastic adaptors such as those supplied with the product, we recommend using ~8-12mm wide PTFE tape, wrapped around the male thread before tightening.
- 3 Ensure that the system is correctly connected. The 'donut' labels around the ports are clearly marked with inlet and outlet symbols and function in both English and French language. Ports marked as outlet means fluid leaves the product and must be connected to the process inlet.
- 4 Check all joints are tight and leak free.
- Where this product is incorporated into other equipment, it is the responsibility of the assembler to ensure safety



Installation, Operation & Service Manual **G04/K04**

Rev 1 - 05/04/23

4.3 BACKFILLING

- In situations where the chiller is situated physically lower than the application being cooled, fluid will apply pressure to the water circuit of the product.
- 2 The weakest seal is normally the tank lid, and this is typically where fluid will escape the unit.
- 3 Ideally, the product will be located higher or level with the product waterline. If this is not possible, a non-return solenoid valve kit can be installed as an optional standard assembly.
- 4 Please raise any questions with the sales team on sales@app-therm.com.

4.4 **VOLTAGE SELECTION**



Different chillers will have different requirement for electrical supply. The specification and requirement can be found on the back plate of the chiller. Please ensure that the voltage supply and number of phases are correct on the site and your power supply can handle the stated current draw.

4.4 FILLING PROCEDURES

- a) Check all application valves are open, including solenoid valves and variable position valves.
 - i) The product will require an open water circuit to pump into.
 - ii) Any obstructions can increase the time, or entirely prevent the bleeding of air from the system.
- b) Remove the cap from the tank lid on the top of the product.
 - i) Fill the tank to just underneath the bottom of the filling port.
- c) Switch the unit on
 - i) Immediately after toggling this switch (assuming the power cord is connected and switched on), the pump will start running.
 - ii) Keep the unit running until the water level in the tank drops and cuts out on the level switch interlock.
 - iii) Turn the unit off at the main switch.
 - iv) Fill the tank again to lift the level switch.
 - v) Start the product



Caution: Do not run the pump dry. Do not deadhead the pump.

d) Repeat step c) until the chiller does not cut out on the level switch.

e) Now with the unit running:

- i) Leave the cap off the tank for >30mins to allow air to escape, or very loosely screwed on to prevent water splashing out of the tank if the unit has a 'flow through' design.
- f) Check the application and tubing for signs of leaks whilst the chiller is running.
 - i) Replace the tank lid fully when satisfied the system is full and bled of air and no abnormal sound can be heard from the pump

Installation, Operation & Service Manual **G04/K04**

Rev 1 - 05/04/23

4.5 DRAINING A SYSTEM

- a) Isolate the unit. Have a suitable bucket or drain on-hand.
 - i) Remove the tank lid to allow air into the system.
 - ii) Disconnect hoses individually. Be aware hoses as well as the unit are filled with fluid.
 - iii) Use tank drain tap (If fitted)
 - iv) Consider using the red transport plugs to block product ports to give yourself time to empty hoses before continuing to empty the product.
- b) Local rules affect where fluid can be disposed of. Ensure hazardous products do not enter the water course and are reclaimed from the unit for professional disposal.



Installation, Operation & Service Manual **G04/K04**

Rev 1 - 05/04/23

5.0 OPERATION

KO4/GO4 series chillers are fitted with a high performance KR1 PID controller, which can control temperature to within 0.1½C of set point. In addition, there is a high and low temperature warning via the LED on the display of the controller, which is triggered if the temperature deviates more than 10½C from the set point.

102C	from the set point.							
	DISPLAY CONTENTS DURING NORMAL OPERATION							
1	Physical navigation buttons, up, down, return and enter.							
2	8888.8 is the actual read value on input sensor.							
3	888.8 is the setpoint value.							
4	Rectangles bottom left 1-4 display when output is active.							
5	MAN LED shows in manual mode (fixed output value).							
6	°C or °F shows units as settable in the 'inP' group.							
7	AL LED appears when output is beyond a set alarm point.							
	ACCESS TO SETTINGS							
8	Push the return button for more than 5 seconds. The upper display will show PASS while the							
	lower display will show 0.							
9	Using up and down buttons set the programmed password – full access is granted by entering							
	'40'. ATC are not responsible for damage either to the chiller or the connected equipment as a							
	result of changing parameters without ATC's oversight.							
10	During parameter modification the instrument continues to perform process control. In certain							
	conditions, when a configuration change can produce a significant change to the process, it is							
	advisable to temporarily stop the controller from controlling during the programming procedure							
	(control outputs will be OFF).							
	A password equal to 2000 + the programmed value (i.e. 2000 + 40 = 2040). The control will							
	restart automatically when the configuration procedure will be manually closed.							
11	Push the return button. If the password is correct the display will show the acronym of the first							
	parameter group 'inP'. Push button for more than 5 seconds, the instrument will come back to							
	the "standard display".							
12	The configuration parameters are collected in various groups. Every group defines all parameters							
42	related with a specific function (control, alarms, output functions).							
13	Push return button for more than 5 seconds, the instrument will come back to the "standard							
	display". For specific settings and guidance, review the controller datasheets provided by ATC. If							
	you are not in receipt of these, please contact ATC using the information in the header of this document.							
	GENERAL NAVIGATION							
14	Return button; A short press allows to exit from the current parameter group and select a new							
	parameter group. A long press allows you to close the configuration parameter procedure (the							
	instrument will come back to the "standard display").							
15	Enter button; When the upper display is showing a group and the lower display is blank, this key							
	allows to enter in the selected group. When the upper display is showing a parameter and the							
	lower display is showing its value, this key allows to store the selected value for the current							
	parameter and access the next parameter within the same group.							
16	Up button; Allows to increase the value of the selected parameter.							
17	Down button; Allows to decrease the value of the selected parameter.							
18	Pushing both Return and Enter buttons moves back to the previous group. Press return first to							
	start. The selection of the group is cyclic (on a carousel), so it is possible to move back around to							
	the group you require.							

Installation, Operation & Service Manual **G04/K04**

Rev 1 - 05/04/23


5.1 CHANGING SET POINT

	ACHIEVING DESIRED SET POINT
1	Push and hold the enter button until it shows "SP 1"
2	Use the up and down arrows to set the temperature
3	Press enter button again to set the temperature
4	The set value now will be activated, and the chiller will start

5.2 REGULATING PRESSURE/FLOW

POSITIVE DISPLACEMENT PUMPS & DISCRETE PRESSURE RELIEF VALVE

This arrangement comprises a positive displacement pump (most commonly a rotary vane type) with a spring-loaded pressure relief valve to provide better overpressure setting control with minimal flow losses compared to fixed orifice bypasses/reliefs. This annex describes ATC's default settings and how to adjust the system.




It is important to understand the basic principle that all else being equal, higher flow results in a higher demand for pressure. The motor generates the power required to turn the pump head and create that pressure. The more restrictive a water circuit is, the higher the pressure required to maintain flowrate. Positive displacement pumps are designed to generate high pressure and are mechanically tight - their RPM dictates flowrate. See image above.

- **Pump discharge** fluid leaves the pump head and enters the gauge-tee assembly. It passes the temperature sensor at the tee, and when the pressure relief is inactive, fluid heads towards no. 2.
- **Outlet to process** fluid heads out of the product to the application. Whatever restrictions lie downstream, fluid leaving this point has not passed through the pressure relief valve.
- **Pressure gauge connection** connection to a pressure gauge on the front panel of the product. This displays the pressure in the water circuit, at the pump outlet, not the application which will be lower.
- 4 Pressure Relief Valve (PRV) spring tension in the body of the PRV determines whether fluid travels through point 2 or passes through to point 7.

Installation, Operation & Service Manual **G04/K04**

Rev 1 - 05/04/23

5	PRV adjustment knob – rotate clockwise to increase maximum delivery pressure. Rotate anti-clockwise to limit maximum delivery pressure.
6	PRV adjustment locknut – release this locknut to adjust no.5. Ensure it is tightened once adjustments are made. If it is not, vibration may cause the adjustment knob to move on itself.
7	Bypass flow outlet – where the pressure requirement to overcome restrictions downstream of no.2 rises to be higher than the setting at no.5, the spring inside no.4 will compress and allow liquid to start bleeding through to no.7. It is important to understand that the nature of the spring means there can be no black and white point for pressure relief setting – the spring will slowly compress and bleed flow until all flow passes through no.7. When fully bypassing, all flow stays inside the chiller to protect the application.



Installation, Operation & Service Manual **G04/K04**

Rev 1 - 05/04/23

6.0 BASIC TROUBLESHOOTING

	AIRCOOLED CHILLER TROUBLESHOOTING								
#	SYMPTOM	CAUSES							
1	Compressor or pump motor not running.	Level switch may have been tripped – Fill the tank up again to resolve the issue High pressure switch may have been tripped:							
	Large condenser fan is running, and control illuminated.	 a) Check fan spins freely and that condenser is clear of debris. b) Assess whether the ambient temperature is too high. c) Assess whether ventilation/air circulation is poor. d) Assess whether surrounding equipment is generating hot air and 							
		feeding it into the chiller air intakes. If fitted, a low-pressure switch may have been tripped – this indicates refrigerant loss. A multimeter is required to test this. If running at a high setpoint, or if operating conditions have changed and raised the fluid temperature, it is possible the compressor has turned itself on high temperature internal protection.							
2	Noisy Operation	Pump motor or fan motor bearings may have failed. Pump head vane or seal may have failed. In a 3-phase system, phase rotation may be incorrect.							
3	High fluid pressure/low flow rate	Fluid lines may have become fouled through solid debris or biological growth. Ensure biocides or pre-mixed biocidal process fluids are being used. Ensure hoses are not transparent or translucent – UV light entering prompts growth. General fluid leak can lead to flow rate failure. Excess fittings, too many sharp direction changes, too small a diameter hose for flowrate or hoses too long at the wrong diameter.							
4	Fluid collecting or leaking	Leaks are impossible to diagnose remotely, but usually you will see a pool of fluid on a lower surface. Do not confuse this with condensation from cold parts of the fridge. Review the chemical compatibility of your fluid with the chiller. Contact ATC if you are unsure.							
5	Poor Cooling Performance	Check if the fans are running. Clean the condenser coil. Excess thermal load applied. Assess whether the ambient temperature is too high. Assess whether ventilation/air circulation is poor. Assess whether surrounding equipment is generating hot air and feeding it into the chiller air intakes.							
	Please contact ATC if	you could not find the information you were looking for or still have an issue with the unit							



Installation, Operation & Service Manual **G04/K04**

Rev 1 - 05/04/23

7.0 END USER MAINTANANCE

Caution: Failure to carry out service at the specified intervals may permanently													
damage your equipment.													
Print this sheet out and display close to the chiller to maximize the visibility of maintenance requirements.													
Weekly									4				
Check Fluid Level – Top up as required													
Monthly	J	F	'	М	Α	M	J	J	Α	S	0	N	D
Check the condenser is free from dust or accumulation of debris													
. "				Y1	,	Y2	Y3	, [Y4		Y5	,	Y6
Annually				11		1 2	13	•	14		10		10
Drain process flui replace with fresh													
Check for fluid le													
throughout chiller and application													
Clear any debris from inside the chiller													
A vacuum cleaner is recommended for cleaning out the condenser, while soft cloths and IPA are recommended for cleaning metallic surfaces. If any spillages have occurred, best practice is to allow the water to evaporate off and wipe up remaining glycol residue with a cloth. Always clean with power supply isolated.													
Caution: N	Caution: Never blow out the condenser with compressed air												
Caution: If						amage	d, con	tact A	TC who	will b	e able	to sup	ply



Installation, Operation & Service Manual **G04/K04**

Rev 1 - 05/04/23

8.0 DECLARATION AND APPROVALS



Operating Manual; Declarations & Approvals

Annex J-5

Date 03/FEB/2021 Author(s) MJH Page 1/1 Revision 01

CONFLICT MINERALS COMPLIANCE STATEMENT

Applied Thermal Control (ATC) adheres to and embraces the ethical values that support our everyday activities. As an expression of these principles and ethical values, ATC adheres to the principle of responsible sourcing of components containing precious and non-precious metals and metal salts in compliance with applicable laws and regulations.

The metals considered are Tantalum (Ta), Tungsten (W), Tin (Sn) and Gold (Au). ATC actively sources components from suppliers known to be reputable and could demonstrate compliance upon request with the Conflict Minerals acts and guidelines.

ATC uses Gold and Tin in electrical components, on PCBs and in rotating machinery, as governed by technical requirements of products. These metals could potentially originate from conflict mineral sites. As many of our suppliers do not purchase these metals direct from smelters, both they and ATC must rely heavily on information that will be provided by their suppliers to determine the source and chain of the metals in those products.

ATC is committed to working with its customers and supply chain to meet the customer's specification and requirements with regards to traceability, sourcing requirements and restrictions. ATC commits that, to the best of our knowledge, our suppliers are complying with the conflict minerals act as stated in their documentation. These statements are reviewed, and updates obtained as required.

Robert Poniatowski, CEO Signed in Barrow-upon-Soar, UK, date 03/FEB/2021



Operating Manual; Declarations & Approvals

Annex J-7

DOCUMENT DETAILS

Date 03/FEB/2021

Author(s) MJH

Page 1/1

Revision 01

WHAT IS THE REACH REGULATION 1907/2006?

REACH is a regulation of the European Union, adopted to improve the protection of human health and the environment from the risks that can be posed by chemicals, while enhancing the competitiveness of the EU chemicals industry. REACH places the burden of proof on companies. To comply with the regulation, companies must identify and manage the risks linked to the substances they manufacture and market in the EU. They have to demonstrate to ECHA how the substance can be safely used, and they must communicate the risk management measures to the users. If the risks cannot be managed, authorities can restrict the use of substances in different ways. In the long run, the most hazardous substances should be substituted with less dangerous ones. REACH stands for Registration, Evaluation, Authorization and Restriction of Chemicals. It entered into force on 1/JUN/2007.

REACH 'ARTICLE' COMPLIANCE CONSIDERATIONS

REACH ANNEX XVII COMPLIANCE

Substances under Annex XVII are restricted either in full (not to be used at all) or for specific uses (can be used in some uses but cannot be used in identified uses).

Applied Thermal Control has contacted all our suppliers and to the best of our knowledge, none of the articles that we sell intentionally contain any of the Annex XVII substances currently on the candidate list in concentrations of >0.1% by weight.

REACH ANNEX XIV COMPLIANCE

Substances under Annex XIV require authorization to use in the EU after sunset date, require communication to downstream recipients when over threshold (0.1% w/w at article level) and require notification to ECHA when SVHC over threshold and imported over 1000kg annually and use not already registered.

Applied Thermal Control has contacted all our suppliers and to the best of our knowledge, none of the articles that we sell intentionally contain any of the Annex XVII substances currently on the candidate list in concentrations of >0.1% by weight.

SVHC LIST COMPLIANCE

Substances of Very High Concern (SVHC) require communication to downstream recipients when over threshold (0.1% w/w at the article level), notification to the European Chemicals Agency (ECHA) when SVHC over threshold and when imported over 1000kg annually and use not already registered.

Applied Thermal Control has contacted all our suppliers and to the best of our knowledge, none of the articles that we sell intentionally contain any of the Annex XVII substances currently on the candidate list in concentrations of >0.1% by weight.

DECLARATION

Robert Poniatowski, CEO Signed in Barrow-upon-Soar, UK, date 03/FEB/2021



Operating Manual; Declarations & Approvals

Annex J-8

DOCUMENT DETAILS

Date 03/FEB/2021 Author(s) MJH Page 1 / 1 Revision 01

WHAT IS THE POPs REGULATION 2019/1021?

POPs stands for persistent organic pollutants. In Europe, the global Stockholm Convention is implemented through POPs legislation. POPs are organic substances that persist in the environment, accumulate in living organisms and pose a risk to our health and the environment. They can be transported by air, water or migratory species across international borders, reaching regions where they have never been produced or used. International risk management is necessary as no region can manage the risks posed by these substances alone.

The European Parliament (and Council) issued regulation 2019/1021 on 20/JUN/2019, and further amended (regulation 2020/784) on 8/APR/2020.

POPs LISTED UNDER INITIAL REGULATION 2019/1021

Pesticides:

Aldrin, Chlordane, DDT, Dieldrin, Endrin, Heptachlor, Hexachlorobenzene, Mirex, Toxaphene.

Industrial Chemicals;

Hexachlorobenzene, Polychlorinated Biphenyls (PCBs).

Industrial Chemical Byproducts;

Hexachlorobenzene byproducts;

Polychlorinated dibenzo-p-dioxins and polychlorinated dibenzofurans (PCDD/PCDF), and PCBs.

POPs LISTED UNDER AMENDMENT 2020/784

Perfluorooctanoic acid (PFOA), its salts and PFOA-related compounds.

POPS COMPLIANCE STATEMENT

We certify that to the best of our knowledge, based upon up-to-date information from our suppliers, all products supplied by Applied Thermal Control are fully POPs compliant in accordance with regulations and amendments above mentioned.

DECLARATION

Robert Poniatowski, CEO Signed in Barrow-upon-Soar, UK, date 03/FEB/2021



Operating Manual; Declarations & Approvals

Annex J-10

DOCUMENT DETAILS

Date 03/FEB/2021 Author(s) MJH Page 1 / 1 Revision 02

WHAT IS THE ROHS DIRECTIVE?

The RoHS Directive places restrictions on the use of certain hazardous substances in electrical and electronic equipment (EEE). RoHS compliance has been required for many years, however in 2014 it became a mandatory requirement under CE Marking. ATC products do not clearly fall within any of the existing categories of equipment, but as of 23/JUL/2019, all EEE not covered falls within scope of the directive. In contrast to RoHS 1, RoHS 2 is a CE marking Directive, and requires, for finished EEE, the use of the CE mark on the product to show compliance. The responsibility for affixing the CE mark resides with the manufacturer.

RoHS 1 2002/95/EC

Adopted in February 2003 by the EU and taking effect on 1/JUL/2006, RoHS 1 restricted the use of 6 hazardous materials;

- 1) Lead (Pb)
- 2) Mercury (Hg)
- 3) Cadmium (Cd)
- 4) Hexavalent Chromium (Cr6+)
- 5) Polybrominated Biphenyls (PBB)
- 6) Polybrominated Diphenyl Ether (PBDE)

We certify that to the best of our knowledge, based upon up-to-date information from our suppliers, all products supplied by Applied Thermal Control are fully RoHS 1 compliant.

RoHS 2 2011/65/EU

Adopted in July 2011 by the EU and taking effect on 2/JAN/2013, RoHS 2 expands the scope of RoHS 1 by adding new categories. RoHS 2 compliance is required to CE mark the product. Compliance with RoHS 2 is mandatory from 22/JUL/2019.

We certify that to the best of our knowledge, based upon up-to-date information from our suppliers, all products supplied by Applied Thermal Control are fully RoHS 2 compliant.

RoHS 3 2015/863/EU

Adopted in 2015 by the EU and taking effect from 22/JUL/2019, RoHS 3 adds four additional substances to RoHS 1's list.

- 1) Bis(2-Ethylhexyl) phthalate (DEHP): < 1000 ppm
- 2) Benzyl butyl phthalate (BBP): < 1000 ppm
- 3) Dibutyl phthalate (DBP): < 1000 ppm
- 4) Di-isobutyl phthalate (DIBP): < 1000 ppm

We certify that to the best of our knowledge, based upon up-to-date information from our suppliers, all products supplied by Applied Thermal Control are fully RoHS 3 compliant.

DECLARATION

Robert Poniatowski, CEO Signed in Barrow-upon-Soar, UK, date 03/FEB/2021



Operating Manual; Declarations & Approvals

Annex J-2

DOCUMENT DETAILS

Date 6/APR/2022 Author(s) MJH Page 1 / 1 Revision 1

EU DECLARATION OF CONFORMITY

Document layout; Governed by Machinery Directive 2006/42/EC, Annex II.

REGISTERED BUSINESS ADDRESS

Applied Thermal Control Ltd, 39 Hayhill Industrial Estate, Barrow-upon-Soar, Loughborough, LE12 8LD, UK.

AUTHORISATION TO COMPILE THE TECHNICAL FILE

Mitchell Howard, address as above

DESCRIPTION & IDENTIFICATION OF MACHINERY

	DESCRIPTION & IDENTIFICATION OF MACHINER I
Generic denomination;	K-Series
Function;	Recirculating chiller
Model;	All with 'K' prefix.
Type;	Air-cooled or water-cooled vapour compression-based.
Serial number;	
Commercial name;	As above.

NOTIFIED BODY

Not applicable

QUALITY ASSURANCE SYSTEM

QMS International Ltd, Muspole Court, Muspole Street, Norwich, NR3 1DJ, UK. ASCB Registered; 201409-2

DECLARATION

Applied Thermal Control declares that the machinery described above fulfils all the relevant provisions of the directives and standards below.

Directive	Harmonised Standards applied
Machinery Directive 2006/42/EC (inclusive Low Voltage Directive 2014/35/EU)	EN ISO 12100:2010 (MD) BS EN 61010-1:2010+A1:2019 (LVD)
EMC Directive 2014/30/EU	IEC 61000-6-2:2005 IEC 61000-6-4:2006 +A1:2011
RoHS Directive 2011/65/EU (RoHS 2) RoHS Directive (EU) 2015/863 (RoHS 3)	EN IEC 63000:2018
Pressure Equipment Directive (2014/68/EC)	Out of Scope. Sound Engineering Practice (SEP) applied.

PERSON EMPOWERED TO DRAW UP DECLARATION

Robert Poniatowski, CEO

Signed in Barrow-upon-Soar, UK, date 6/APR/2022



Operating Manual; Declarations & Approvals

Annex J-20

DOCUMENT DETAILS

Date 20/DEC/2022 Author(s) WSE Page 1 / 1 Revision 1

UKCA

UKCA DECLARATION OF CONFORMITY (DoC)

Demand The Product Safety and Metrology etc. (Amendment etc.) (EU created by; Exit) Regulations 2019

REGISTERED BUSINESS ADDRESS

Applied Thermal Control Ltd, 39 Hayhill Industrial Estate, Barrow-upon-Soar, Loughborough, LE12 8LD, UK.

AUTHORISATION TO COMPILE THE TECHNICAL FILE

Mitchell Howard, Applied Thermal Control Ltd, 39 Hayhill Industrial Estate, Barrow-upon-Soar, Loughborough, LE12 8LD, UK.

DESCRIPTION & IDENTIFICATION OF MACHINERY

K-Series
Water Heat Exchanger
All with 'K' prefix.
Water cooled heat exchanger.
As above.

NOTIFIED BODY

Not applicable

QUALITY ASSURANCE SYSTEM

QMS International Ltd, Muspole Court, Muspole Street, Norwich, NR3 1DJ, United Kingdom. ASCB Registered; 201409-2

DECLARATION

The manufacturer declares that the machinery described above is in conformity with the relevant statutory requirements applicable to the specific product. The manufacturer takes full responsibility for the product's compliance.

- Supply of Machinery (Safety) Regulations 2008
- Electromagnetic Compatibility Regulations 2016
- Electrical Equipment (Safety) Regulations 2016
- The Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012

PERSON EMPOWERED TO DRAW UP DECLARATION

Robert Poniatowski, CEO Signed in Barrow-upon-Soar, UK, date 20/DEC/2022

Issue 6.4, January 2023 Page 1 of 4



Conforming to Directive 1907/2006/EC

SECTION 1: IDENTIFICATION	NOF THE SUBSTANCE/MIXTURE AND O	F THE COMPANY/UNDERTAKING

1.1. Product Name Hexid A4

1.2. Manufacturer Applied Thermal Control Limited

39 Hayhill Industrial Estate, Barrow upon Soar, Leicestershire, LE12 8LD. United Kingdom.

www.app-therm.com

1.3. Telephone Number +44(0)1530 839998
 1.4. Email sales@app-therm.com
 1.5. Emergency Telephone Number +44(0)1530 839998
 1.6. Intended/Recommended Use Heat Transfer Fluid

SECTION 2: HAZARDS IDENTIFICATION

2.1. Classification of the substance or mixture

The product is not classified as dangerous according to Regulation (EC) No. 1272/2008. This mixture is not classified as dangerous according to Directive 1999/45/EC.

SECTION 3: COMPOSITION/INFORMATION ON INGREDIENTS

3.1. Chemical Nature Water (CAS 7732-18-5), not classified.

Propylene glycol (CAS 57-55-6) (REACH 01-2119456809-23)

(EINECS 200-338-0) not classified.

Fluorescein (trace) and biocide (trace) not classified.

3.2. Food Grade

SECTION 4: FIRST AID MEASURES

	General advise	No special precautions required. Treat symptomatically.
4.1.	Eye Contact	Rinse thoroughly with plenty of water, also under the eyelids. Remove contact lenses
		after a few minutes and continue rinsing. If symptoms persist, call a physician.
4.2.	Skin Contact	Wash off immediately with plenty of water. If skin irritation persists, call a physician.
4.3.	Inhalation	Remove to fresh air. If symptoms persist, call a physician.
4.4.	Ingestion	Rinse mouth with water. Never give anything by mouth to an unconscious person. If
		symptoms persist, call a physician.

SECTION 5: FIREFIGHTING MEASURES

5.1. Extinguishing media

Use extinguishing measures that are appropriate to local circumstances and the surrounding environment. Water spray, foam, dry powder or CO2. Alcohol-resistant foam

5.2 Unsuitable extinguishing Media

High volume water jet. Do not use a solid water stream as it may scatter and spread fire.

5.3 Specific hazards during firefighting

In fire conditions, toxic decomposition products may be formed (see also section 10). In combustion, emits fumes, smoke, carbon dioxide (CO2) and carbon monoxide (CO). Heating will cause a pressure rise - with severe risk of bursting and explosion, Violent steam generation or eruption may occur upon application of direct water to hot liquids.

5.4. Advice for firefighters

In the event of fire, wear self-contained breathing apparatus. Wear personal protective equipment. Promptly isolate the scene by removing all persons from the vicinity of the incident if there is a fire. Keep containers cool by spraying with water if exposed to fire. Collect contaminated fire extinguishing water separately. This must not be discharged into drains. Burning fluids may be extinguished by dilution with water

Issue 6.4, January 2023 Page 2 of 4

Conforming to Directive 1907/2006/EC

SECTION 6: ACCIDENTAL RELEASE MEASURES

Personal precautions 6.1.

Use personal protective equipment. Avoid contact with skin and eyes. Keep unnecessary and unprotected personnel from entering the area.

6.2. Precaution to protect the environment

Do not flush into surface water or sanitary sewer system. Avoid subsoil penetration.

6.3 Clean-up procedures

> Contain the spillage, soak up with non-combustible absorbent material, (e.g. sand, earth, diatomaceous earth, vermiculite) and transfer to a container for disposal according to local / national regulations (see section 13). Keep in suitable, closed containers for disposal. Dike the area of spill to prevent spreading and pump liquid to salvage tank. Treat recovered material as described in section 13 Disposal considerations.

SECTION 7: HANDLING AND STORAGE

Precautions for safe handling

Keep container tightly closed. Handle in accordance with good industrial hygiene and safety practice. Spills of these organic materials on hot fibrous insulations may lead to lowering of the auto-ignition temperatures possibly resulting in spontaneous combustion.

7.2. Conditions for safe storage

Keep only in the original container.

SECTION 8: EXPOSURE CONTROLS/PERSONAL PROTECTION

8.1. Control parameters

Material

Component: Propane-1.2-diol CAS-No. 57-55-6

Other Occupational Exposure Limit Values EH40 WEL, Time Weighted Average (TWA):, Total vapour and particulates.150 ppm, 474 mg/m3

EH40 WEL, Time Weighted Average (TWA):, Particulate.10 mg/m3

ELV (IE), Time Weighted Average (TWA):, Total vapour and particulates.150 ppm, 470 mg/m3

ELV (IE), Time Weighted Average (TWA):, Particulate.10 mg/m3

8.2. Exposure controls/Appropriate engineering controls

> Local exhaust. If this product contains ingredients with exposure limits, use process enclosures, local

exhaust ventilation or other engineering controls to keep worker exposure below any

recommended or statutory limits.

Personal protective equipment

Respiratory protection Suitable respiratory protective device Combination filter: A-P2

Filter Type Combined particulates and organic vapour type

Hand protection Category short time exposure Break through time> 10 min

Protective index Class 1 When prolonged exposure is expected: Break through time> 120 min Protective index Class 4 Observe the information of the glove manufacturers on permeability. Protective gloves should be chosen according to Workplace Safety Assessment.

Gloves recommended according to EN 374 (protection against chemicals). Chemical resistant gloves made of butyl rubber or nitrile rubber category III

according to EN 374.

SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES

9.1	Appearance at 20°C	Fluorescent green clear liquid
9.2	Odour	Almost odourless
9.3	Flash point	Boils without flashing
9.4	Ignition temperature	Not Available
9.5	Flammability Limit	Not Available
9.6	Oxidizing Properties	Not Available
9.7	Auto flammability	450°C
9.8	Density at 25°C	~1.036g/cm ³
9.9	pH (as is)	7
9.10	Boiling point	102°C
9.7	Auto flammability	450°C
9.8	Solubility in water	Miscible
9.9	Freezing point	-21°C
9.10	Specific Heat Capacity	3.78kJ/kg °K
9.11	Viscosity, Kinetic, at 25°C	3.51mPa.s

Issue 6.4, January 2023 Page 3 of 4

Conforming to Directive 1907/2006/EC

SECTION 10: STABILITY AND REACTIVITY

10.1. Reactivity

Stable under recommended storage conditions. No dangerous reaction known under conditions of normal use.

10.2. Chemical stability

No decomposition if stored and applied as directed. Stable under recommended storage conditions. Hygroscopic.

10.3. Hazardous reactions

Hazardous polymerisation does not occur.

10.4. Conditions to avoid

Generation of gas from decomposition causes pressure in closed systems. Keep away from direct sunlight. Avoid high temperatures. Avoid temperatures exceeding the decomposition temperature. Avoid UV light.

10.5. Materials to avoid

Strong acids, Strong bases, Strong oxidizing agents.

10.6. Hazardous decomposition products

Aldehydes, Alcohols, Ether, Organic acids.

SECTION 11: TOXICOLOGICAL INFORMATION

11.1. Toxicity Oral

LD50: > 20000 mg/kg (rat) This product can present a small hazard if large quantities are swallowed.

11.2. Inhalation

LC50 : 6.15 mg/l (rat; 4 h; vapour) At ambient temperature the exposure to vapours is minimal due to a low volatility rate. Inhalation may cause irritation to the nose, throat, upper respiratory tract and lungs. No deaths occurred

11.3. Dermal

LD50: > 20000 mg/kg (rabbit) Prolonged skin contact is unlikely to result in absorption of harmful amounts. Skin irritation by prolonged exposure is unlikely. Repeated contact may cause flaking and softening of skin.

11.4. Eyes

Slight irritation is possible. Direct contact with eyes may cause temporary irritation. Corneal injury is unlikely.

11.5. Sensitisation

Patch test on human volunteers did not demonstrate sensitisation properties.

11.6. CMR Carcinogenicity

Animal testing did not show any carcinogenic effects. Information given is based on data obtained from similar substances.

11.7. Mutagenicity

No data available.

11.8. Reproductive toxicity

No data available.

11.9. Specific Target Organ Toxicity

Single exposure no data available. Repeated exposure no data available.

11.10. Other toxic properties

Repeated dose toxicity. In rare cases, repeated excessive exposure to propylene glycol may cause central nervous system effects. Aspiration hazard Due to its physical properties, the substance does probably not pose any aspiration hazard.

11.11. Other relevant toxicity information

Handle in accordance with good industrial hygiene and safety practice.

11.12. Experience with human exposure

Health injuries are not known or expected under normal use.

Conforming to Directive 1907/2006/EC

SECTION 12: ECOLOGICAL INFORMATION

12.1. Acute toxicity

Fish - LC50: 40613 mg/l (Oncorhynchus mykiss; 96 h) (static test)

Daphnia and other aquatic invertebrates - LC50: 18340 mg/l (Ceriodaphnia Dubia (water flea); 48 h) (static test)

Algae - ErC50 : 19000 mg/l (Pseudokirchneriella subcapitata (green algae); 96 h) (Growth inhibition)

Bacteria - NOEC: > 20000 mg/l (Pseudomonas putida; 18 h)Chronic toxicity

Aquatic invertebrates - NOEC: 13020 mg/l (Ceriodaphnia Dubia (water flea); 7 d) (semi-static test)

12.2. Persistence and degradability

Biodegradability 81 % (anaerobic; Exposure Time: 28 d)(OECD 301 F)

Readily biodegradable 96 % (anaerobic; Exposure Time: 64 d)(OECD 306.)

12.3. Bioaccumulative potential

BCF - 0.09 estimated Low bioaccumulative potential

12.4. Mobility

Estimated Koc < 1, indicating very high soil mobility.

12.5. PBT and vPvB assessment

Not a PBT or vPvB substance or mixture

12.6. Other adverse effects

Do not flush into surface water or sanitary sewer system. Avoid subsoil penetration. This substance is not in Annex I of Regulation (EC) 2037/2000 on substances that deplete the ozone layer.

SECTION 13: DISPOSAL CONSIDERATION

13.1. Waste treatment methods

Disposal together with normal waste is not allowed. Special disposal required according to local regulations. Do not let product enter drains. Contact waste disposal services.

13.2. Contaminated packaging

Empty contaminated packaging thoroughly. They can be recycled after thorough and proper cleaning. Packaging that cannot be cleaned are to be disposed of in the same manner as the product.

13.3. European Waste Catalogue Number

No waste code according to the European Waste Catalogue can be assigned for this product, as the intended use dictates the assignment. The waste code is established in consultation with the regional waste disposer.

SECTION 14: TRANSPORT INFORMATION

Not dangerous goods for ADR, RID, IMDG and IATA.

14.1. EEC Regulations

UNNO None Class None Packing Group None

Road & Rail Transport (ADR & RID) None IMDG Not Applicable ICOA None

14.2. Export commodity code

39074000

Classification - Polycarbonates.

14.3. Weight and dimensions

5Kg per 5 litre container. 19x14x29cm.

14.4. Manufactured in the United Kingdom

SECTION 15: REGULATORY INFORMATION

15.1 Classification Not classified as hazardous to users.

15.2. CAS No. 5755615.3. Risk or Safety phrases None15.4. Labelling None

SECTION 16: OTHER INFORMATION

Key literature references and sources for data taken from supplier information and data from the "Database of registered substances" of the European Chemicals Agency (ECHA) were used to create this safety data sheet. Other information - The information provided in this Safety Data Sheet is correct to our knowledge at the date of its revision. The information given only describes the products with regard to safety arrangements and is not to be considered as a warranty or quality specification and does not constitute a legal relationship.

The information contained in this Safety Data Sheet relates only to the specific material designated and may not be valid for such material used in combination with any other material or in any process, unless specified in the text.

W HARP
INTERNATIONAL

HARP[®] 407C

Revision: 1

Revision date: May 2020

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier

Product name	R407C
REACH registration number	See section 3: Composition/information on ingredients
CAS No.	See section 3: Composition/information on ingredients
EC No.	See section 3: Composition/information on ingredients

1.2 Relevant identified uses of the substance or mixture and uses advised against

Product use	Refrigerant Industrial uses: Uses of substances as such or in preparations at industrial sites Professional uses: Public domain (administration, education, entertainment, services, craftsmen)
Restricted use	Consumer uses: Private households (= general public = consumers)
Description	Gas

1.3 Details of the supplier of the safety data sheet

	210 Details of the supplier of the surety acts sheet	
Company	Harp International Limited	
Address	Gellihirion Industrial Estate	
	Pontypridd	
	Rhondda Cynon Taff	
	CF37 5SX	
	UK	
Web	www.harpintl.com	
Telephone	+44 (0) 1443 842 255	
Fax	+44 (0) 1443 841 805	
Email	harp@harpintl.com	
Email of competent person	safety@harpintl.com	

1.4 Emergency telephone number

Emergency telephone number	+44 (0) 1270 502 891
	24 hours

SECTION 2: Hazards identification

2.1 Classification of the substance or mixture

Classification – EC 1272/2008	Compressed gas: H280
-------------------------------	----------------------

2.2 Label elements

2.2 Label elements	
Hazard pictograms	
Signal word	Warning
Hazard statement	H280 – Contains gas under pressure; may explode if heated
Precautionary statement	P410+P403 – Protect from sunlight. Store in a well-ventilated place.

2.3 Other hazards

according to Regulation (EO) 2013



HARP[®] 407C

Revision: 1

Revision date: May 2020

SECTION 3: Composition/information on ingredients

3.1 Substances

EC 1272/2008

Chemical name	CAS No.	EC No.	REACH registration	Concentration	Classification
			number	(%w/w)	
1,1,1,2-Tetrafluoroethane (R134a)	811-97-2	212-377-0	01-2119459374-33	ca. 52	Compressed gas: H280
Pentafluoroethane (R125)	354-33-6	206-557-8	01-2119485636-25	ca. 25	Compressed gas: H280
Difluoromethane (R32)	75-10-5	200-839-4	01-2119471312-47	ca. 23	Flam. Gas 1: H220 Compressed gas: H280

SECTION 4: First aid measures

4.1 Description of first aid measures

Inhalation	Move the exposed person to fresh air
Eye contact Rinse immediately with plenty of water	
Skin contact	Frostbite: treat as thermal burns
Ingestion	Ingestion is not considered a potential route of exposure

4.2 Most important symptoms and effects, both acute and delayed

Inhalation	Symptoms may include loss of mobility/consciousness. Victim may not be
	aware of asphyxiation. Respiratory arrest.
Eye contact	Contact with liquefied gas can cause damage due to evaporative cooling
Skin contact	Contact with liquefied gas can cause damage due to evaporative cooling
Ingestion	Ingestion is not considered a potential route of exposure

4.3 Indication of any immediate medical attention and special treatment needed

Inhalation	If you feel unwell, seek medical advice
Eye contact	Seek medical attention if irritation or symptoms persist
Skin contact	Thaw frosted parts with lukewarm water. Do not rub affected area. Get immediate medical advice/attention.
	ininiediate inedical advice/attention.
Ingestion	Ingestion is not considered a potential route of exposure

SECTION 5: Firefighting measures

5.1 Extinguishing media

	This	product	is	not	flamı	mable	in	air	under	ambient	conditi	ions	of
	temp	perature	and	pres	sure.	Use	extin	guis	hing m	edia app	ropriate	to	the
	surro	ounding fi	re c	ondit	ions.								
E O Consider the second section of the const	.1	L -4		•	_								

5.2 Special hazards arising from the substance or mixture

At high temperature, toxic and/or corrosive fumes may be produced by
thermal decomposition (gaseous hydrogen fluoride (HF), carbon oxides).

5.3 Advice for firefighters

Wear self-contained breathing apparatus and protective clothing. Heat may cause the containers to explode. Keep fire exposed containers cool by
spraying with water. Fire exposed containers may vent contents through pressure relief devices. In case of fire nearby, remove exposed containers.

SAFETY DATA SHEET

according to Regulation (EU) 2015/830



HARP[®] 407C

Revision: 1

Revision date: May 2020

SECTION 6: Accidental release measures

6 1 Dorconal	Inrocautions	protoctivo		224	emergency procedures
6.1 Personal	i precautions.	. protective	eaulbment	and	emergency procedures

Ensure adequate ventilation of the working area. Avoid contact with skin
and eyes. Evacuate personnel to a safe area. Wear self-contained breathing
apparatus and protective clothing. Vapours are heavier than air. Prevent
from entering sewers, basements or workpits. Do not enter confined spaces
where gas may have accumulated.

6.2 Environmental precautions

Prevent further leakage or spillage if safe to do so.

6.3 Methods and material for containment and clean up

Allow to evaporate. Provide adequate ventilation.

6.4 Reference to other sections

See section 8 Exposure controls / personal protection See section 13 Disposal considerations

SECTION 7: Handling and storage

7.1 Precautions for safe handling

Only experienced and properly instructed persons should handle gases
under pressure. Protect containers from physical damage. Do not drag, roll,
slide or drop. Do not remove or deface labels. Adopt best manual handling
considerations when handling, carrying and dispensing. Secure cylinders in
an upright position at all times. Close valves when not in use and when
empty. Ensure adequate ventilation of the working area. Do not allow
backfeed into the container. Avoid contact with skin and eyes. When using,
do not eat, drink or smoke. Never use direct flame or electrical heating
device to raise the pressure of the container.

7.2 Conditions for safe storage, including any incompatibilities

	Keep containers tightly closed. Keep in a cool, dry, well-ventilated area. Store in correctly labelled containers. Keep away from sources of ignition – no smoking. Store out of direct sunlight.
Suitable packaging	Stainless steel, steel.

7.3 Specific end use(s)

, 10 opening on a doc(s)	
	See section 1.2 Relevant identified uses of the substance or mixture and
	uses advised against for further information.

SECTION 8: Exposure controls/personal protection

8.1 Control parameters – exposure limit values

Component	CAS No.	Value type	Exposure limit values	Source
		(form of exposure)		
1,1,1,2-Tetrafluoroethane	811-97-2	TWA	1000ppm / 4240mg/m ³	EH40
Pentafluoroethane	354-33-6		Not listed in EH40	
Difluoromethane	75-10-5		Not listed in EH40	



Revision: 1

Revision date: May 2020 **8.2 Exposure controls**

HARP[®] 407C

•	
Appropriate engineering controls	Ensure adequate ventilation of the working area. Oxygen detectors should
	be used when asphyxiating gases may be released. Systems under
	pressure should be regularly checked for leaks.
Individual protection measures	Wear protective clothing
Eye/face protection	Approved safety goggles
Skin & body protection	Wear suitable gloves. Wear safety shoes when handling containers.
Respiratory protection	Wear suitable respiratory protection equipment when necessary
Occupational exposure controls	Keep away from food, drink and animal feedstuffs.
Hygiene protection	Good industrial hygiene and safety procedures. Do not eat, drink or smoke
	when using the product.

SECTION 9: Physical and chemical properties

9.1 Information on basic physical and chemical properties

Appearance - Physical state	Gas
Appearance - Form	Liquefied gas
Colour	Colourless
Odour	Slight ethereal
Odour threshold	No data available
рН	Not applicable
Melting point	No data available
Boiling point / range	-43.6°C
Flash point	Not applicable
Evaporation rate	Not applicable
Flammability (solid, gas)	This product is not flammable
Upper explosion limit / Lower	Not applicable
flammability limit	
Vapour pressure	1,190.3 kPa (25°C)
Vapour density	3.03 (air = 1)
Relative density	No data available
Solubility(ies)	
Water solubility	No data available
Partition coefficient:	
n-octanol/water	No data available
Auto-ignition temperature	Not applicable
Decomposition temperature	No data available
Viscosity	
Viscosity, kinematic	No data available
Explosive properties	Not applicable
Oxidising properties	Not applicable

SECTION 10: Stability and reactivity

10.1 Reactivity

,	Stable under normal conditions
10.2 Chemical stability	
	Stable under normal conditions. The gaseous product in the presence of air can form, under certain conditions of temperature and pressure, a flammable mixture.

acco. a....B to 110Batation (207 202



HARP[®] 407C

Revision: 1

Revision date: May 2020

10.3 Possibility of hazardous reactions

	No data is available on this product	
10.4 Conditions to avoid		
	Keep away from heat and sources of ignition. Avoid contact with flames and red hot metallic surfaces.	
10.5 Incompatible materials		
	Alkaline hydroxides, alkaline earth metals, strong oxidizing agents, finely	
	divided metals.	
10.6 Hazardous decomposition products		
	Under normal conditions of storage and use, hazardous decomposition products should not be produced.	
	At high temperature, thermal decomposition can give rise to toxic and corrosive products.	

SECTION 11: Toxicological information

11.1 Information on toxicological effects

1111 mornidation on toxicological circuit	
Acute toxicity	As with other volatile aliphatic halogenated compounds, through vapour accumulation and/or inhalation of large quantities, the product can cause loss of consciousness and cardiac disorders aggravated by stress and lack of oxygen. Risk of mortality.
Skin corrosion/irritation	Ejection of liquefied gas: frostbite possible
Serious eye damage/irritation	Ejection of liquefied gas: frostbite possible
Respiratory or skin sensitisation	No data available
Germ cell mutagenicity	No data available
Carcinogenicity	No data available
Reproductive toxicity	No data available
STOT single exposure	No data available
STOT repeated exposure	No data available
Aspiration hazard	No data available
Repeated or prolonged exposure	No data available

SECTION 12: Ecological information

12.1 Toxicity

	No data available	
12.2 Persistence and degradability		
	Not applicable to gases and gas mixtures	
12.3 Bioaccumulative potential		
	Expected to biodegrade and not expected to persist for long periods in an	
	aquatic environment	
12.4 Mobility in soil		
	Unlikely to cause ground or water pollution due to its high volatility	
12.5 Results of PBT and vPvB assessment		
	Not classified as PBT or vPvB	
12.6 Other adverse effects		
	Contains fluorinated greenhouse gases. When discharged in large quantities may contribute to the greenhouse effect. Global warming potential: 1774	

W HARP
INTERNATIONAL

HARP[®] 407C

Revision: 1

Revision date: May 2020

SECTION 13: Disposal considerations

13.1 Waste treatment methods

Dispose of in accordance with all local and national regulations. Avoid
discharges to atmosphere. Refer to manufacturer/supplier for information
on recovery/recycling. Dispose of container via supplier only.
EWC code: 14 06 01* Chlorofluorocarbons, HCFC, HFC

SECTION 14: Transport information

Hazard pictograms



14.1 UN number

	UN 3340

14.2 UN proper shipping name

REFRIGERANT GAS R407C

14.3 Transport hazard class(es)

14.5 Transport nazara ciass(cs)	
ADR/RID	
Class	2
Labels	2.2
Hazard No. (ADR)	20
Tunnel category	(C/E)
Emergency action code	2TE
IMDG	
Class	2.2
EmS No.	F-C, S-V
IATA	
Class	2.2
Packing instruction	200

14.4 Packing group

- ···· · · · · · · · · · · · · · · · ·		

14.5 Environmental hazards

Environmental hazards	Not applicable
Marine pollutant	Not classified as a marine pollutant

14.6 Special precautions for user

Avoid transport on vehicles where the load space is not separated from
the driver's compartment. Ensure container valves are closed, not leaking
and caps in place. Ensure containers are firmly secured. Ensure adequate
air ventilation.

14.7 Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not applicable

WHARPINTERNATIONAL

HARP[®] 407C

Revision: 1

Revision date: May 2020

SECTION 15: Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

Regulations	REGULATION (EC) No 1907/2006 OF THE EUROPEAN PARLIAMENT AND OF					
	THE COUNCIL of 18 December 2006 concerning the Registration,					
	Evaluation, Authorisation and Restriction of Chemicals (REACH),					
	establishing a European Chemicals Agency, amending Directi					
	1999/45/EC and repealing Council Regulation (EEC) No 793/93 and Commission Regulation (EC) No 1488/94 as well as Council Directive					
	76/769/EEC and Commission Directives 91/155/EEC, 93/67/EEC,					
	93/105/EC and 2000/21/EC.					

15.2 Chemical safety assessment

No CSA has been carried out

SECTION 16: Other information

Other information

Text of Hazard Statements in	H280: Contains gas under pressure; may explode if heated.		
Section 3	H220: Extremely flammable gas		
Reference materials	HSE publication EH40/2005 Workplace exposure limits (latest edition)		
Changes from previous versions	-		

Further information

The information supplied in this safety data sheet is designed only as
guidance for the safe use, storage and handling of the product. This
information is correct to the best of our knowledge and belief at the date
of publication however no guarantee is made of its accuracy. This
information relates only to the specific material designated and may not
be valid for such material used in combination with any other materials or
in any other process. Material users should review the information and
recommendations in the specific context of their intended manner of
handling, use, processing and storage, including an assessment of the
appropriateness of the material in the user's end product, if applicable.



Tel: +44(0) 1530 839998

Temperature Control for Industrial, Medical & Scientific Applications

sales@app-therm.com www.app-therm.com

Section 1: Identification of the substance/ mixture and of the company/ undertaking

1.1 Product identifier

Product Name CoolFlow EG

1.2 Relevant identified uses of the substance or mixture and uses advised against

Identified uses Industrial grade secondary refrigerant with antifreeze and inhibitor functions

designed for use in process cooling, refrigeration and HVAC systems.

Uses advised against This product is not recommended for any industrial, professional or consumer use

other than the identified uses above.

1.3 Details of the supplier of the Safety Data Sheet

Supplier Applied Thermal Control Ltd

39 Hayhill Industrial Estate

Barrow Upon Soar leicestershire LE12 8LD

+44 (0) 01530 839998 sales@app-therm.com

1.4 Opening Hours Monday-Thursday - 08:00 - 17:00 Friday - 08:00 - 14:00

Section 2: Hazards identification

2.1 Classification of the substance or mixture

Classification - Regulation (EC) No. 1272/2008 (CLP)

Physical and chemical hazards
Human health
Environment

Not classified as a physical or chemical hazard
Acute Tox. 4 - H302, STOT RE 2 - H373
Not classified as an environmental hazard

2.2 Label elements

EC No. N/A

Labelling - Regulation (EC) No. 1272/2008 (CLP)

Pictograms





Signal Word - Warning



Applied Thermal Control Ltd 39 Hayhill Industrial Estate Barrow Upon Soar

Leicestershire, LE12 8LD

Tel: +44(0) 1530 839998

Temperature Control for Industrial, Medical & Scientific Applications

sales@app-therm.com www.app-therm.com

Hazard statements

H302 - Harmful if swallowed

H373 - May cause damage to organs - Kidneys - through prolonged or repeated exposure if swallowed.

Precuationary statements

P260 - Do not breathe dust/fumes/gas/mist/vapours/spray

P264 - Wash hands thoroughly after handling

P270 - Do not eat, drink or smoke when using this product

P301+P312 - IF SWALLOWED: Call a POISON CENTRE or doctor/physician if you feel unwell

P301+P330+P331 - IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.

2.3 Other hazards

This product does not meet the PBT/vPvB criteria of REACH, annex XIII.

Section 3: Composition/information on ingredients

3.2 Mixtures

Component - Monoethylene glycol (ethane-1, 2-diol)

 Concentration
 80-95%

 EC No.
 203-473-3

 CAS No.
 107-21-1

Reach registration No. 01-2119456816-28

Classification - Regulation (EC) No. 1272/2008 (CLP)

Acute Tox. 4 - H302 STOT RE 2 - H373

Component - Ethanol

 Concentration
 1-5%

 EC No.
 200-578-6

 CAS No.
 64-17-5

Reach registration No. 01-211945719-43

Classification - Regulation (EC) No. 1272/2008 (CLP)

Flam. Liq. 2 - H225 Eye Irrit. 2 - H319

Component - Sodium nitrite

 Concentration
 <0.25%</td>

 EC No.
 231-555-9

 CAS No.
 7632-00-0

Reach registration No. 01-2119471836-27

Classification - Regulation (EC) No. 1272/2008 (CLP)

Ox. Sol. 3 - H272 Acute Tox. 3 - H301 Eye Irrit. 2 - H319 Aquatic Acute 1 - H400



Tel: +44(0) 1530 839998

Temperature Control for Industrial, Medical & Scientific Applications

sales@app-therm.com www.app-therm.com

Section 4: First aid procedures

4.1 Description of first aid procedures

General Information When safe to do so remove the victim from the source of exposure giving

consideration as to whether this may cause further discomfort to the victim.

Inhalation Move the affected person to fresh air at once. Keep warm in a position

comfortable for breathing. If breathing becomes difficult, properly trained personnel may assist the victim by supplying oxygen to ease breathing.

Get medical attention if any discomfort continues.

Ingestion Do NOT induce vomitting and seek medical attention immediately. Move victim to

fresh air and keep warm and at rest in a position comfortable for breathing. Give victim approximately 250 mL of water however, do not give victim anything to drink if not fully conscious. If medical advice is delayed and an adult has consumed several ounces of this chemical, give approximately 100 mL of hard liquor (for

children give 2 mL per kilogram of body weight).

Skin Contact Remove contaminated clothing immediately and wash skin with soap and water.

Get medical attention if any discomfort continues.

Eye Contact Rinse immediately with plenty of water. Remove any contact lenses and open

eyelids wide apart. Continue to rinse for at least 15 minutes.

Get medical attention if any discomfort continues.

4.2 Most important symptoms and effects, both acute and delayed

General Information The following symptoms are listed in case of exposure to the 100% neat product.

Inhalation Ingestion

Skin Contact

Inhalation of vapours may cause mild irritation of the upper respiratory tract. Initial symptoms may include an upset stomach, nausea, vomiting and diarrhoea. Symptoms may progress to hyperventilation, metabolic acidosis, cardiovascular

dysfunction and acute kidney failure depending on the extent of poisoning. Prolonged and repeated contact may cause mild irritation of the skin.

Eye Contact Direct eye contact may cause reddening of the eyes.

4.3 Indication of immediate medical needs or special treatment

If several ounces (> 50 mL) of this product have been ingested, early administration of ethanol may help to counteract the toxic side effects such as metabolic acidosis, cardiovascular dysfunction and in severe cases kidney failure. Consider haemodialysis or peritoneal dialysis and thiamine 100 mg plus pyridoxine 50 mg intravenously every 6 hours. If ethanol is used a therapeutically effective blood concentration in the range of 100-150 mg/dL may be achieved by a rapid loading dose followed by a continuous intravenous fusion.

Respiratory symptoms such as pulmonary edema, may be delayed. Victims receiving significant exposure should be kept under observation for 24-48 hours for signs of respiratory distress. In the case of severe poisoning, respiratory support with mechanical ventilation and oxygenation of the patient.

Notes for the doctor

No specific recommendations in addition to the suggestions in Sections 4.1 and 4.3. Treat symptomatically.



Tel: +44(0) 1530 839998

Temperature Control for Industrial, Medical & Scientific Applications

sales@app-therm.com www.app-therm.com

Section 5: Firefighting measures

5.1 Extinguishing media

Extinguish with alcohol-resistant foam, carbon dioxide (CO2), dry chemicals, sand and dolomite or water fog.

5.2 Special hazards arising from the substance or mixture

Specific Hazards When heated and in the case of a fire, harmful vapours/gases (such as carbon

monoxide and carbon dioxide) may be formed.

Unusual fire and explosion hazards Exposure to extreme heat may cause product containers to explode.

5.3 Advice for firefighting

Protective actions during firefighting Move containers away from fire area if this can be done without risk.

Keep people away, isolate the fire and deny unnecessary entry.

Use water fog to keep fire-exposed containers cool and disperse vapours. Runoff water should be prevented from entering sewers and watercourses.

Specialist protective equipment for Wear positive-pressure self-contained breathing apparatus (SCBA) and full

Fire fighters protective clothing.

Section 6: Procedure for unwanted emissions

6.1 Personal precautions, protective equipment and emergency procedures

Personal precautions

Avoid flames, sparks, heat and smoking.

In the case of inadequate ventilation, use respiratory protection.

Protective Equipment

Wear protective clothing as described in Section 8 of this Safety Data Sheet.

Emergency Procedures

Stop leak/release if possible to do so without risk.

Extinguish all ignition sources if safe to do so.

Warn everybody of potential danger and evacuate if necessary.

6.2 Environmental precautions

Do not discharge into drains, water courses or onto the ground.

Spillages or uncontrolled discharges into watercourses must be reported immediately to the Environmental Agency or other appropriate regulatory body.

6.3 Methods and materials for containment and clean up

Absorb spillage with inert, damp, non-combustible material, then flush the contaminated area with water. Containers with collected spillage should be appropriately labelled with the correct contents and hazard labels. Collect and place in suitable waste disposal containers and seal securely. For waste disposal, see Section 13.

6.4 Reference to other sections

Wear protective clothing as described in Section 8 of this Safety Data Sheet. Collect and dispose of spillage as indicated in Section 13.



Tel: +44(0) 1530 839998

Temperature Control for Industrial, Medical & Scientific Applications

sales@app-therm.com www.app-therm.com

Section 7: Handling and storage

7.1 Precautions for safe handling

Avoid spilling and contact with the skin and the eyes as well as direct inhalation of sprays and mists. Provide good ventilation.

Do not eat, drink or smoke in work areas and wash hands after handling this product.

7.2 Conditions for safe storage including any compatibilities

Store in tightly-closed, original containers.

Keep separate from food, feedstuffs, fertilisers and other sensitive material.

Do not store near heat sources or expose to high temperatures.

Keep away from heat, sparks and open flame.

7.3 Specific end use(es)

The identified uses for this product are detailed in Section 1.2.

Section 8: Exposure controls / Personal protection

8.1 Control parameters

Name	STD	TWA-8 Hrs	STEL-15 Min
Monoethylene glycol (ethane-1, 2-diol)	WEL	52 mg/m ³	104 mg/m ³

DNEL

Industry, Inhalation - Long term local effects: 35mg/m³ Industry, Dermal - Long term systemic effects: 106mg/m³ Consumer, Inhalation - Long term local effects: 7mg/m³ Consumer, Dermal - Long term systemic effects: 7mg/m³

PNEC

Fresh water: 10 mg/L Marine water: 1mg/L STP: 199.5 mg/L

Sediment fresh water: 20.9 mg/kg

Soil: 1.53 mg/kg

NameSTDTWA-8 HrsSTEL-15 MinEthanolWEL1920 mg/m³Not available

DNEL

Industry, Inhalation - Short term local effects: 1900 mg/m³ Industry, Dermal - Long term systemic effects: 343 mg/kg/day Industry, Inhalation - Long term systemic effects: 950 mg/m³ Consumer, Inhalation - Short term local effects: 950 mg/m³

Consumer, Dermal - Long term systemic effects: 206 mg/kg/day Consumer, Inhalation - Long term systemic effects: 114 mg/m³ Consumer, Oral - Long term systemic effects: 87 mg/kg/day

PNEC

Fresh water: 0.96 mg/L

Marine water: 0.79 mg/L



Tel: +44(0) 1530 839998

Temperature Control for Industrial, Medical & Scientific Applications

sales@app-therm.com www.app-therm.com

Intermittent release: 2.75 mg/L

STP: 580 mg/L

Sediment fresh water: 3.6 mg/kg Sediment marine water: 2.9 mg/kg

Soil: 0.63 mg/kg

Name STD TWA-8 Hrs STEL-15 Min Sodium nitrite WEL Not available Not available

DNEL

Industry, Inhalation - Long term systemic effects: 2 mg/m³ Industry, Inhalation - Short term systemic effects: 2 mg/m³

PNEC

Fresh water: 0.0054 mg/L Marine water: 0.00616 mg/L Intermittent release: 0.0054 mg/L Sediment (fresh water): 0.0195 mg/kg Sediment (marine water): 0.0223 mg/kg

STP: 21 mg/L

8.2 Exposure controls



Technical procedures

Other Protection

Engineering measures Methods to prevent or control exposure are preferred. Provide adequate ventilation to minimise the risk of inhalation of sprays and mists.

Handle in accordance with good industrial hygiene and safety practices. Hygiene measures

Wash hands after handling this product and at the end of each work shift. Routinely wash work clothing and personal protective equipment to remove

possible contaminants.

Respiratory equipment Where risk assessment shows air-purifying respirators are appropriate use a full-

> face respirator with multi-purpose combination (US) or type ABEK (EN 14387) respirator cartridges as a backup to engineering controls. If the respirator is the sole means of protection, use a full-face supplied air respirator. Use respirators and components tested and approved under appropriate government standards

such as NIOSH (US) or CEN (EU).

PVC/butyl rubber/neoprene gloves are recommended. Hand protection

Eye protection Wear approved chemical goggles or face shield. Skin Protection Wear rubber apron or protective clothing in case of contact.

Wear suitable protective clothing/footwear as protection against splashing or

contamination.



Tel: +44(0) 1530 839998

Temperature Control for Industrial, Medical & Scientific Applications

sales@app-therm.com www.app-therm.com

Thermal Hazards No specific measures required.

Environmental Exposure Controls Product not classified as an environmental hazard - no specific environmental

exposure controls required.

Section 9: Physical and chemical properties

9.1 Information on basic physical and chemical properties

Blue liquid **Appearance**

Odour Non-pungent but characteristic aroma

Not applicable Odour Threshold

7.5 - 10.5 depending on inhibitor formulation

-40°C Melting point / Pour point >180°C Initial boiling point Flash point >64°C

Evaporation Rate No test data available

Product is not classified as flammable Flammability

Flammability / explosion limits Upper limit: 13% Lower limit: 3%

Vapour pressure No test data available Vapour density (air = 1) No test data available

Relative density of the mixture 1.04 - 1.2 Solubility Soluble in water Partition coefficient: n-octanol / water No test data available

Auto-ignition temperature >390°C

Decompostion temperature No test data available Viscosity See product data sheet

Explosive properties Not applicable - product is not classified as an explosive Oxidising properties Not applicable - product is not classified as an oxidising agent

9.2 Other information

Not determined.

Section 10: Stability and reactivity

10.1 Reactivity

There are no known reactivity hazards associated with this product.

10.2 Chemical stability

Stable at normal ambient temperatures and when used as recommended.

Product is hygroscopic and will absorb water by contact with the moisture in the air.

10.3 Possibility of hazardous reactions

There are no known hazardous reactions associated with this product.

10.4 Conditions to avoid

Avoid temperatures >180°C for prolonged periods of time, flames and sources of ignition.

10.5 Incompatible materials

Strong acids, strong alkalis and strong oxidising agents.

10.6 Hazardous decomposition products



Tel: +44(0) 1530 839998

Temperature Control for Industrial, Medical & Scientific Applications

sales@app-therm.com www.app-therm.com

No known hazardous decomposition products. Potentially hazardous products released due to fire are listed in Section 5.2 of this Safety Data Sheet.

Section 11: Toxicological information

11.1 Information on toxicological effects

Acute toxicity

This product has not been tested as a whole for acute toxicity health effects. For this reason, the acute toxicity values for the main components of this mixture are listed below.

Acute toxicity values for monoethylene glycol:

LD50, oral, rat: 4700 mg/kg bw

LD50, dermal, rabbit: >10,600 mg/kg bw LD50, dermal, mouse: >3500 mg/kg bw LC50, inhalation (vapours), rat: >2.5 mg/L

Oral toxicity of monoethylene glycol is expected to be moderate in humans despite animal testing showing a lower degree of toxicity. The estimated lethal dose in humans of NEAT MONOETHYLENE GLYCOL is expected to be approximately 100mL.

Acute toxicity values for ethanol: LD50, oral, rat: 10,470 mg/kg bw LD50, dermal rabbit: 17,100 mg/kg bw

Skin corrosion/ irritation
Serious eye damage/ irritation
Skin irritation is not expected when this product is used/handled correctly.
Eye irritation is not expected when this product is used/handled correctly.

Respiratory/ skin sensitisation

Germ cell mutagenicity
Carcinogenicity

Carcinogenicity Product is not expected to be carcinogenic.

Percentage to demand the representation of the product is not expected to demand the representation.

Reproductive toxicity Product is not expected to damage the reproductive system or harm a developing

Product is not expected to be mutagenic.

Product not classified as a skin/respiratory sensitiser.

fetus.

Evaluation of CMR properties

STOT-single exposure

No test data available. No test data available.

STOT-repeated exposure

Aspiration hazard

No test data available. No test data available.

General information

See Section 4.2 of this Safety Data Sheet.

Inhalation Inhalation of vapours may cause mild irritation of the upper respiratory tract.

Ingestion Initial symptoms may include an upset stomach, nausea, vomiting and diarrhoea.

Symptoms may progress to hyperventilation, metabolic acidosis, cardiovascular dysfunction and acute kidney failure depending on the extent of poisoning.

Skin contact Prolonged and repeated contact may cause mild irritation of the skin.

Eye contact Direct eye contact may cause reddening of the eyes.



Tel: +44(0) 1530 839998

sales@app-therm.com

www.app-therm.com

Temperature Control for Industrial, Medical & Scientific Applications

Section 12: Ecological information

Ecotoxicity

The product is not classified as hazardous to the environment.

12.1 Toxicity

LC50, 96 hours, fish: >100 mg/L - not classified as harmful to fish >100 mg/L - not classified as harmful to daphnia EC50, 96 hours, aquatic plants: >100 mg/L - not classified as harmful to aquatic plants

12.2 Persistence and degradability

This product is readily biodegradable (90% over 10 days).

12.3 Bioaccumulative potential

Will not bio-accumulate.

Partition coefficient - not determined.

12.4 Mobility in soil

Product is mobile in soil as it is water soluble.

12.5 Results of PBT and vPvB assessment

This product does not meet the PBT/vPvB criteria of REACH, annex XIII.

12.6 Other adverse effects

Not determined.

Section 13: Advice on disposal

General information

Waste to be treated as controlled waste. Disposal to licensed waste disposal site in accordance with Local Waste Disposal Authority.

Disposal methods

Dispose of waste and residues in accordance with local authority and/or local sewage treatment plant requirements.

Section 14: Transport information

14.1 UN number

Product not hazardous for transport - no information required.

14.2 UN proper shipping name

Product not hazardous for transport - no information required

14.3 Transport hazard class(es)

Product not hazardous for transport - no information required.

Transport labels

Product not hazardous for transport - no information required.

14.4 Packing group

Product not hazardous for transport - no information required.



Tel: +44(0) 1530 839998

sales@app-therm.com

Temperature Control for Industrial, Medical & Scientific Applications

www.app-therm.com

14.5 Environmental hazards

Product not classed as an environmentally hazardous substance or marine pollutant.

14.6 Special precautions for user

Product not hazardous for transport - no information required.

14.7 Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code

Product not hazardous for transport - no information required.

Section 15: Regulatory information

15.1 Safety, health and environmental regulations / legislation for the substance or mixture

EU legislation

Regulation (EC) No 1907/2006 of the European Parliament and of the Council of 18 December 2006 concerning the Registration, Evaluation Authorisation and Restriction of Chemicals (REACH) (as amended). Regulation (EC) No 1272/2008 of the European Parliament and of the Council of 16 December 2008 on classification, labelling and packaging of substances and mixtures (as amended).

Guidance notes

CHIP for everyone HSG228. Approved Classification and Labelling Guide (Sixth edition) L131. Safety Data Sheets for substances and preparations.

15.2 Chemical safety assessment

No chemical safety assessment for this mixture has been carried out.

Section 16: Other information

Issued by Thermal Exchange Ltd Revision Date 21/03/2018 Approved by Graham Wade Revision Comments Review in line with CLP regulation

Hazard statements in full

The following hazard statements are the hazard statements 'in full' for the components of this mixture. They are not the hazard statements associated with the final classification of this product.

H302 - Harmful if swallowed

H373 - May cause damage to organs - Kidneys - through prolonged or repeated exposure if swallowed

H225 - Highly flammable liquid and vapour

H319 - Causes serious eye irritation

H272 - May intensify fire; oxidiser

H301 - Toxic if swallowed

H400 - Very toxic to aquatic life

Further classification and composition comments

No further classification or composition comments required.



Tel: +44(0) 1530 839998

Temperature Control for Industrial, Medical & Scientific Applications

sales@app-therm.com www.app-therm.com

(i) Indication of changes

Safety Data Sheet updated to comply with the new requirements as set out in Regulation (EC) No. 1272/2008 (CLP).

(ii) Abbreviations and acronyms

bw: bodyweight

CAS No: Chemical Abstracts Service number

CLP: Classification Labelling and Packaging Regulation

DNEL: Derived No-Effect Level EC: European Commission

EC No: European Chemical number: EINECS, ELINCS or NLP

ECHA: European Chemicals Agency

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

LC50: Lethal Concentration, 50%

LD50: Median Lethal Dose

PBT: Persistent, Bioaccumalative & Toxic PNEC: Predicted No Effect Concentration

REACH: Registration, Evaluation, Authorisation & restrictions of Chemicals

SDS: Safety Data Sheet

vPvB: Very Persistent and Very Bioaccumalative

WEL: Workplace Exposure Limit

(iii) Training advice

Product should only be handled by trained operators.

(iv) Additional information

The above information describes exclusively the safety requirements of the product and is based on our present-day knowledge. The information is intended to give advice about the safe handling of the product named in this Safety Data Sheet, for storage, processing, transport and disposal. The information cannot be transferred to other products. In the case of mixing the product with products or in the case of processing, the information on this Safety Data Sheet is not necessarily valid for the new made-up material.